

Improve Data Quality!

Live high resolution data viewing immediately alerts you to any problems (e.g. poor electrode contact), so they can be fixed as they arise- resulting in higher quality recordings!

Easily connect to:

- A standalone PC (Using MA-720 & WINDAQ-PRO)
- Directly to an exiting ADC (Such as Motion Capture)
- Simultaneously Connect to Both (Using free splitter cable)




Analog to Digital Data Conversion Hardware (MA-720)



- Connects MA300 EMG System to a PC for stand-alone EMG!
- Easy to Connect and Use- USB Connection
- High Resolution Capability- (250 kHz at 16-bits)


Real Time Viewing Software (WINDAQ-LITE)



- Included with the MA-720 HARDWARE
- Real Time Heads-Up Display- View, live streaming data.
- Display Up to 16 EMG channels.
- Multitasking operation while still displaying live data!

Analog to Digital Conversion Hardware.....MA-720 (Includes WINDAQ-LITE)

Advanced Data Recording Software (WINDAQ-PRO)



- Record and display 32 data channels, live and in real time.
- Sample your data at optimal speeds! Up to 250,000/sec.
- Export, import & translate data files (text, ASCII, MATLAB...)
- Label channels and add data comments during/after collection.
- Includes many data handling and analysis features such as Automatic Time and Date Stamping, Over-Sampling, Frequency Analysis, Digital Filtering, X-Y Plotting, and Programmable Gain.

Advanced Data Recording Software.....WINDAQ-PRO

EASY INTEGRATION- Works with VICON, MAC, Qualysis, CodaMotion, AMTI, Bertec, Kistler, and more!

Motion Lab Systems, Inc.
www.motion-labs.com

Baton Rouge, Louisiana, USA

Email: sales@motion-labs.com

Ph: (225) 272-7364 Fax: (225) 272-7336