

Flexibility! With 3 different channel types available in the MA300 EMG System you can easily connect a wide range of equipment.

Goniometers

We've made it easy for our customers to get synchronized goniometer and torsionmeter data. Ideal for quick, simple, and accurate measurement of joint movement in multiple planes -

Features

- Extremely robust, lightweight and flexible
- Undetectable- comfortably worn under clothing
- Does not hinder the actual movement of the joint

Goniometer Interface Box (Specify Connector).....GIB-000
Twin Axis Electro Goniometers (Five size options).....SG Series
Single Axis Gonio/Torsionmeters (Two sizes options).....Q Series
Mini Single Axis Electro Goniometer (e.g. fingers).....F Series
Goniometer Cables at 1500mm (60in).....A-1500
Goniometer Extension Cable (GIB to MA300).....A-1SGQ-002


Event Switches

Easily Detect Event Contact.

This collection of three separate items can be used together for easy connection to Motion Lab Systems EMG Systems or used independently for custom designs.

Features

- Accurately detect contact
- Small 15mm on/off sensor
- Saves Money- reusable low cost design
- Comfortable- 0.5mm ultra thin switch
- Easy placement- flexible 100mm tail
- Easy to use in custom applications



A. Event Switch-25mm dia. w/ 100mm Tail.....MA-153
B. Event Switch Extension Cords (pack of 8).....MA-136
C. Event Switch Adaptor Cable Pair.....Part Codes Vary
Double Sided Sticky Tape (1 roll).....000-100-TPE

www.motion-labs.com

Baton Rouge, Louisiana, USA

Email: sales@motion-labs.com

Ph: (225) 272-7364 Fax: (225) 272-7336