

## Subject Vests and Belts

Neoprene Subject Vest and Belts are custom made by Motion Lab Systems to comfortably and securely hold the MA300 EMG Systems and other equipment in place.

### Features

- Secure- Large Velcro placement pad to securely hold the EMG System backpack in place
- Organized- Included Velcro cable holders help group and organize cables
- Washable and Reusable
- Adjustable- Velcro tab for easy size adjustment
- Safe- Latex free



#### Neoprene Subject Vest

- A. Child Size Vest.....MA-142
- B. Teen Size Vest.....MA-144
- C. Adult Medium Size Vest.....MA-145
- D. Adult Large Size Vest.....MA-143

#### Neoprene Subject Belt

- E. Child Size Belt.....MA-140
- F. Adult Size Belt.....MA-141



## Attachment Straps

Stretchable Straps made of soft fabric with built-in stretch featuring easy size adjustments to securely hold Pre-Amplifiers or other equipment in place without constricting the muscle and still allowing the skin to breath.

### Features

- Washable and Reusable- Neoprene or Velcro Brands
- Adjustable- Velcro tab for easy size adjustment
- Safe- Latex free
- Variety- A range of sizes and can also be cut for custom fits.
- Two inches wide.

	<b>NEOPRENE</b>	<b>VELSTRETCH</b>
--	-----------------	-------------------

25cm (10in)	.....001-010-NEO	.....001-010-VEL
50cm (20in)	.....001-020-NEO	.....001-020-VEL

**FYI-** Two vests are included with the Advanced EMG System,

**Motion Lab Systems, Inc.**  
[www.motion-labs.com](http://www.motion-labs.com)

Baton Rouge, Louisiana, USA

Email: [sales@motion-labs.com](mailto:sales@motion-labs.com)

Ph: (225) 272-7364 Fax: (225) 272-7336

# Disposable Wraps and Tapes

## Coban Wrap

A latex-free disposable wrap tape. This wrap is self adhesive without sticking to skin! Great for securing preamplifiers! A 3M Brand Product. Two inches in width, tan color.



## Double Sided Sticky Tape

Clear disposable double sided sticky tape is sold as a single roll. Latex-free. Great for securing event switches! 1 inch in width.

## Retro-Reflective Tape

This 1" wide single sided adhesive Retro-Reflective Tape by 3M Scotchlite is motion camera compatible and easy to cut. (Sold by the metre or whole rolls available.)



<i>Latex Free Coban Sports Wrap - (One Roll).....</i>	<i>MA-395-CBN</i>
<i>Double Sided Sticky Tape, 25.4mm (1 Inch) Wide - (One Roll).....</i>	<i>000-100-TPE</i>
<i>Scotchlite Retro-Reflective Tape (1 Inch) Wide - (36 inches).....</i>	<i>01-300-TPE</i>


# Electrodes

**Fine Wire Needle Electrodes** are sterilized, disposable, preassembled, stainless steel fine wires. The bi-polar/paired needles have a 200mm long wire for EMG applications and work great with Motion Lab Systems brand discrete pre-amplifiers. (10 pack of individually wrapped assemblies.)




**Mini Adhesive Electrodes** are small single use gel electrodes that adhere directly to the skin surface. Available with a connector or snaps. Great for use with small or hard to reach muscles.

**Reusable Ground Electrode** is a durable reusable electrode with a 1 inch surface area. Connects to the EMG System backpack for an additional system ground.



**Snap Leads** allow you to easily connect snap-on gel electrodes to screw on discrete preamplifiers. (Preamps with built-in snap also available.)

<i>Paired Fine Wire Needle Electrode- Pack of 10 (30mm 27 Gauge).....</i>	<i>000-318-30</i>
<i>Paired Fine Wire Needle Electrode- Pack of 10 (50mm 25 Gauge).....</i>	<i>000-318-50</i>
<i>Cabled Disposable Mini Adhesive Gel Electrodes (or Ground) (3 per pack).....</i>	<i>000-390</i>
<i>Single Ag/AgCl Gel Electrodes- w/ snap connector - (box of 50).....</i>	<i>GELS-101</i>
<i>Dual Ag/AgCl Gel Electrodes- 20mm spacing w/ snap (box of 120).....</i>	<i>GELS-N120.1</i>
<i>Reusable Ground Electrode.....</i>	<i>000-391</i>
<i>Snap leads- Set of Two, each 100mm (4in) long.....</i>	<i>EME-24-SBK</i>



# Bead Markers Set

Includes 16 reusable bead identifiers for easy cable recognition. Snaps on and off.

*Bead Markers Set (16 bead identifiers, 8 colors).....MA-321*

GOT AN IDEA? Call today, new product ideas are always welcome!

